

600 HCN Cold Probe Default Parameters

- **Tuning for Ubiquitin**
 - $^1\text{H} \sim 55$
 - $^{13}\text{C} \sim 65$
 - $^{15}\text{N} \sim 95$
- **Offsets**
 - tof ~ 215.2 (~ 4.77 ppm at 25 deg C)
 - dof (^{13}C) ~ 11509 (~ 174 ppm)
 - dof2 (^{15}N) ~ 1264 (~ 119 ppm)
- **Pulse widths**
 - ^1H : pw = 7.2 at tpwr = 58
 - ^{13}C : pwC = 12.5 at pwClvl = 59
 - ^{15}N : pwN = 31.0 at pwNlvl = 60
- **Deuterium decoupling using channel 4 (assumes garp1 or waltz with pwrD90 of 250 usec)**
 - H2dpwr3D = 44
 - H2dmf3D = 3565
 - H2pw90(CH4) = 146.25 at 49 dB
 - H2seq = garp1
 - H2dres3D = 1
- **Amplifier compressions**
 - compH = 0.97
 - compC = 1.09
 - compN = 0.98
- **Maximum Power Limits (set in ghn_co)**
 - Channel 1: 58dB
 - Channel 2: 60 dB
 - Channel 3: 60 dB
 - Channel 4: 50 dB
- **BioPack Power Limits (set in ghn_co)**
 - BPpwrlimits = 1
 - BPspinlock =
 - ^1H observe = dB
 - ^{13}C decoupling = dB
 - $^1\text{H}/^{13}\text{C}$ spinlock =
 - ^{15}N decoupling = dB
 - ^2D decoupling = dB
- **Gradients**

- gzlvl1 = 18000, gzlvl2 = 18260 (ratio 0.986)
- **Lock level values for lineshape sample**
 - Lock power = 28
 - Lock gain = 32
 - Lock level when shimmed well > 70
- **Gain**
 - Gain values should always be > 18. Values between 20 and 25 are typically good values.